


Electro Stimulation Instrument



Important Notice

- 1) Please read the user manual before any operation.
- 2) Do not put any liquid onto or into the machine. Be careful about the electric connections.
- 3) Avoid operations on pregnant women, or those people with epilepsy, pacemakers or any serious psychological or medical condition, or infections in the area treated
- 4) Please disassemble or assemble the equipment by professional personnel.
Please check if there is any connection inside the machine dropped when some function does not work.

Clean and sterilize the spare parts after each operation.

1, Principle and Effects


By creating electronic pulse to stimulate special points on the body, this kind of function can effectively regulate bio-electric and endocrine system and accelerate metabolism, reaching remarkable effect of weight losing and skin tightening.

The instrument is organically combined with frequency scan circuit and microcomputer system to create stimulant bioelectricity, and accurately holds the essentials of acupuncture of traditional Chinese medicine to make regular deep-seated actuation and massage so as to achieve expedite channels, regularize spirits, break over lymph cyclesystem, remove internal toxin, improve body microcirculation, promote metabolism, excite and redeploy body power against disease, so that the bulk of fat cells can be diminished and the skin can be tightened.

This machine has 10 pairs of electric stimulation pads, which allows operating slimming and skin tightening to the whole body at the same time.

2, Panel and Installation Indication

1) Panel



1. Intensity adjust (L1 to L5)

2. Program select


3. Start/Pause

4. On/OFF

5. Time adjust

6. Intensity adjust (R1 to R5)


7. LCD display



8. Outlet for accessory connection

3, Recommended Operations

1) Turn on the power, and press the button of "4" then the screen will display



2) Connect the accessories as below

① Apply some gel or cream for losing weight on the electro patches and fix them up with the flex on the part of the body.

② Paste the patch on the parts of body as following:


Picture 1: Neck and Shoulder

Function: Stimulate the blood circulation on neck and shoulder which can remove the fat and tire.




Picture 2: Neck and shoulder and arm

Function: Stimulate the blood circulation on neck, shoulder and arm, which can open up to remove fat.




Picture 3: Breast

Function: Stimulate the blood circulation on breast
And stimulate the development of the breast
which can adjust the secretion of the
Hormone.



Picture 4: Belly


Function: Stimulate the blood circulation on belly



and remove the fat and tone up the muscle on chest.


Picture 5 Back

Function: Stimulate the blood circulation on back and tone up the muscle in the back, which can remove the excess fat and alleviate fatigue.



Picture 6: Hip

Function: Stimulate the blood circulation on hip, which remove the excess fat and tone up the muscle of the hip.




Picture 7: Arm

Function: Stimulate the blood circulation on arms and open up channels which can remove the excess fat and tone up the muscle of the arm.



Picture 8 Leg

Function: Stimulate the blood circulation on legs and open up channels which can remove the excess fat and tone up the muscle of the leg.



3) Adjust time and program.

4) Press button of "3" to start

**If any off the intensity adjust "1" and "6" is not turned off before pressing "3", then the machine will alarm with sound, you have to turn it off and then it will start to work.*

4, Packing List

Item Name	Unit	Quantity	Remark
Mainframe	set	1	
Stimulation pad	pair	10	
Pad connector	pce	10	
Elastic band	set	1	
Fuse and power wire	set	1	
Instruction	set	1	

5, Specifications

Voltage	<input type="checkbox"/> 100-110V	<input type="checkbox"/> 220-240V		
Frequency	<input type="checkbox"/> 50HZ	<input type="checkbox"/> 60HZ	Power	35 W
Package	56*42*27 CM		Gross Weight	5 KG